

Continence Formulary (for Adults)

Sussex Community NHS Foundation
Trust

Date of issue: March 2019

Review date: March 2021

South Team

Southfield House

11 Liverpool Gardens

Worthing

01273265855

e-mail sc-tr.continenceteamaaw@nhs.net

This information is intended for use by healthcare professionals within Coastal West Sussex. We have made every effort to check that the information is correct at the time of publication. Coastal West Sussex CCG and Sussex Community Foundation NHS Trust does not accept any responsibility for loss or damage caused by reliance on this information.

Version number 3.1

Continence Formulary

This formulary is a guide for all healthcare professionals that prescribe or request prescriptions for patients that require products due to either bladder or bowel dysfunction.

The intention is that this formulary must be used for all new patients and for current patients when a re-assessment of their needs is completed.

However it is acknowledged that it can be difficult to ensure a uniform approach to the procurement and supply of products across organisations. Whilst prescribing from the formulary is consistent with good clinical practice; it is recognised that there may be occasions when the initiation of a non-formulary product occurs upon discharge from a provider organisation, due to differing procurement arrangements. Whilst there may be good evidence base, safety profile and established practice for the initiation of such product(s), it has to be recognised that a formulary product is the preferred option and therefore review with the view of potentially switching would be expected.

The formulary will be reviewed every two years and at this point any amendments or additions will be considered.

Where ever possible there is a continuity of product use between primary and secondary care, and this has been achieved by joint working on this formulary.

The term "specialist" in this formulary (highlighted in blue, bold and italics) refers to Bladder and Bowel Specialist Nurses. Items with this restricted coding should only be prescribed if recommended following a discussion with the Bladder and Bowel teams.

Version No.	Date	Author(s) of original development or review	Details of document development
1	10/12	Hilary Chiffins – Lead Specialist Continence Nurse	Original Development
2	10/14	Angelique Vaughan – Clinical Nurse Specialist	Full revision
3	6/18	Angelique Vaughan – Clinical Nurse Specialist	Full revision
3.1	6/19	Karen Reynolds – CWS CCG, Senior Pharmaceutical Commissioning Technician	Minor amendments (spelling / format)

This formulary was reviewed by Angelique Vaughan, Lead Specialist Bladder and Bowel Nurse for Sussex Community NHS Foundation Trust in the Coastal West Sussex area, with input from primary care, secondary care and the Medicines Management Team at Coastal West Sussex CCG.

Contents	Page
Urinary Catheters <ul style="list-style-type: none"> • Silicone • Hydrogel Coated Latex • Open ended silicone • Silver alloy coated 	3 - 4
Lubricant gels	4
Catheterisation packs	5
Catheter valves	5
Catheter tube restraining straps	5
Abdominal catheter straps	6
Adhesive catheter clips	6
Catheter maintenance solutions	6
Urinary catheter drainage bags (leg)	6
Leg bag holders	7
Urinary catheter drainage bags(2litre / bed bags)	7
Night bag holders	8
Urinary catheter drainage bags (belly)	8
Miscellaneous supplies	8
Intermittent urinary catheters	9
Urinary sheaths and bags	9
Adhesive free, external collection system	10
Anal Plugs	11
Faecal collector	12
Pads for rectal bleeding	12
Anal irrigation	12
Continence appliances	14
Other useful information/links	15

Urinary Catheters

**Prescribe 3 single units initially, 2 single units on repeat,
Estimated usage 3 single units / 6 monthly**

- Reassess at each catheter change to see if patient is suitable for a TWOC.
- Product choice should be guided by assessment, patient choice and local health economy formulary
- For urethral catheterisation, smallest size that provides adequate draining should be used to avoid bypassing problems.
- Only catheters that have a balloon filled with 10mls of sterile water should be used.
- Standard length urinary catheters must always be used for men. Shorter Female length urinary catheters are suitable for ambulant female patients, but not usually for obese, bed or chair bound female patients. All patients with supra pubic urinary catheters should be re-catheterised with standard length catheters.
- The following urinary catheters are suitable for both urethral and supra pubic catheterisation.
- Long term catheters can stay insitu up to 12 weeks.

Silicone male / standard length (10ml balloon)

Code	Product	Charriere (add size to end of order code)	Ordering
DA3101	Rusch Brilliant Aquaflate	12-16Ch	1 st choice
4833-0512	Great Bear Silicone foley	12-16Ch	2 nd choice
850084-- 000 (add size)	<i>Rusch Brilliant Plus Aquaflate</i> (Cylindrical tip)	<i>12-16</i>	<i>Discuss with specialist</i>

Hydrogel Coated Latex male / standard length (10ml balloon)

Code	Product	Charriere (add size to end of order code)	Ordering
DH3101	Rusch Sympacath Aquaflate	12-16	1 st choice
D2264	Bard Aquamatic Biocath Pre filled with sterile water	12-16	2 nd choice

Silicone female (10ml balloon)

Code	Product	Charriere	Ordering
DA2101	Rusch Brilliant Aquaflate	12-16	1 st choice
AA75	Coloplast Foyzil All Silicone	12-16	2 nd Choice

Hydrogel Coated Latex Female length

Code	Product	Charriere (add size to end of order code)	Ordering
DH210112-	Rusch Sympacath Aquaflate	12-16	1 st choice
D2268	Bard aquamatic Biocath Pre filled with sterile water	12-16	2 nd choice

Open ended silicone Estimated usage 1 single unit / 12 weekly

Suitable for both urethral and supra pubic use and may be a useful alternative for patients who have recurrent bypassing and blocking problems

Code	Product	
<i>AA8C</i>	<i>Coloplast Folsil X-Tra Standard length only</i>	<i>Discuss with specialist</i>

Bard Lubri-sil MALE catheter are suitable choice for Patient who require a comfortable catheter that is hydrogel coated silicone catheter and is useful if patient has latex allergy.

Code	Product	
<i>D1758</i>	<i>Bard Lubri-sil Hydrogel coated silicone</i>	<i>Discuss with Specialist</i>

Bard Lubri-sil FEMALE catheter are suitable choice for Patient who require a comfortable catheter that is hydrogel coated silicone catheter and is useful if patient has latex allergy.

Code	Product	
<i>D1761</i>	<i>Bard Aquafil Lubri-sil Hydrogel coated silicone</i>	<i>Discuss with Specialist</i>

Silver alloy coated Estimated usage 1 single unit / 12 weekly

Urinary catheter reduces bacterial adherence and minimizes biofilm formation.

Code	Product	Charriere (add size to end of order code)
<i>D2365</i>	<i>Bardex I.C. Foley standard length. Male Catheter (28 days only)</i>	<i>12S-22S</i>
<i>D2369</i>	<i>Bardex I.C. Foley Catheter. Female Catheter (28 days only)</i>	<i>12S-16S</i>

Lubricant gels for urinary catheterisations

N.B. The first supra pubic re-catheterisation carried out in the community must be at least 8 weeks post initial, siting carried out in secondary care.

Lubricant with anaesthetic Prescribe 5 as single units / 6 monthly

Lubricants with anaesthetic should be used for both male and female urethral catheterisations except where contra-indicated

Code	Product	Pack size
1160	Optilube active (female)	1 x 6 ml
1161	Optilube active (male)	1 x 11ml
40-006	Instillagel 6ml (female)	10
40-011	Instillagel 11ml (male)	10

Lubricant only **Prescribe 5 as single units / 6 monthly**
Lubricant only is required when carrying our supra-pubic re-catheterisation or where anaesthetic is contra-indicated.

Code	Product	Pack size
1125	Optilube(female)	1x 6ml
1126	Optilube (male)	1 x 11ml

Catheterisation pack **Prescribe 3 packs / 6 monthly**

Two layer system – layer 1 catheter removal kit- with syringe. Layer 2 catheter insertion kit requires syringe for balloon inflation.

Code	Product	Glove size
908410	Cath-it catheterisation packs	S/M
908420	Cath-it catheterisation packs	M/L
908430	Cath- it catheterisation packs	L/XL

Catheter Valves **Estimated usage 1 pack of 5 / 4 weekly**

May be used by patients with long term urethral or supra pubic catheters. Prior to being fitted a risk assessment must be completed and should include the following:

- Does the patient
 - have cognitive awareness to remember to empty their bladder regularly.
 - have the dexterity to use the valve.
 - have an adequate bladder capacity.
- If patient has had surgery recently consultant agreement has been sought.
- Catheter valves must be changed weekly

Code	Product	Pack size
PCV3942	Prosys catheter valve	5
10540A	Great Bear Libra Lever catheter Valve	5

Catheter tube restraining straps **Estimated usage 1 pack of 5 / 4 weekly**

Must be worn by all patients that have a standard urethral urinary catheter insitu and use a leg bag or are bedbound and use a bed bag, to ensure that any traction from the catheter bag is absorbed by the strap. A restraining strap can also help prevent meatal damage. These products are washable and reusable.

Code	Product	Pack size
10644A	Great Bear Fix It (leg strap)	5
408-5536	Proslys catheter strap	short
408-5544		adult
4085551		abdo

Abdominal catheter straps **Estimated usage 1 pack of 5 / 4 monthly**

Abdominal supra pubic catheter tube restraining straps are for use under recommendation of a Bladder and Bowel. They may also be used for patients with high BMI. Straps are washable and reusable.

Code	Product	Pack size
383003	<i>Simpla G Strap (Abdominal strap)</i>	5
10645B	<i>Great Bear Fix It (Abdominal Strap)</i>	5
10647D	<i>Great Bear Fix it (cut to fit length 150 Cm)</i>	5

Adhesive catheter clips **Estimated usage 1 pack of 5 / 4 weekly**

For use in the securing of urethral or suprapubic catheters valves or drainage bags

- Ensure patient reads the removal instruction.
- Patients skin integrity must be careful monitored if this product is used

Code	Product	Pack size
3004	Ugo fix gentle	1 x 5 clips

Catheter maintenance solutions **Estimated usage 1 pack of 10 / PRN**

Only to be used once a clear urinary catheter history is documented and a treatment plan has been written for that named patient. It must be established that the patient may benefit from the planned use of a catheter maintenance solution to extend the life of urinary catheter as per treatment plan, **or** prior to the removal of a urinary catheter to dissolve crystal formation and prevent trauma.

Getcliffe et al (2000) suggest that two sequential instillations of a small volume are more effective than a single administration. BBraun Uro-Tainer Twin, have a twin delivery system allowing for this sequential administration.

Code	Product	Pack size
9746609	Suby G Citric acid 3.23% Twin 2 x 30mls	10
9746625	Solution R Citric acid 6% Twin 2 x 30mls	10
FB99849	Uro -Tainer NaCl 50mls for mechanical flush only (use with Supra pubic catheter change)	10
FB99833	Uro - Tainer NaCl 100mls for mechanical flush only (use with Supra pubic catheter change)	10

Urinary catheter drainage bags **Prescribe 1x10 pack of / 8 weekly**
Leg bags

- Worn either on the thigh or the calf, this will determine the length of the tube.

- There are different types of taps, this must be a considered for patients self-caring.
- Bags need to be changed every 5-7 days. More frequently is not normally necessary and could put the patient at increased risk of infection by breaking the closed system.
- GB Libra conform leg bags have three chambers which minimising audible noises of the urine inside the bag . Weight of urine is more evenly distributes making the Libra more comfortable to wear and less noticeable.

Code	Product	Pack size
Lever action tap (Sterile)		
10114Q	Great Bear Libra Conform leg bag short 500mls	10
10115R	Great Bear Libra Conform leg bag long 500 mls	10
DSL500S	Bard Spirit 500mls short tube 10 cm	10
DSL500L	Bard Spirit 500mls long tube 30 cm	10
Slide Action Tap (Sterile)		
P500S	Prosys Short tube 500mls	10
P500L	Prosys – Long tube 500mls	10

Leg bag holders

Prescribe 1 x pack of 4 / 4 monthly

To hold the urinary catheter leg bag, so reducing the pull on the bladder neck and reducing urethral trauma. Measure the calf or the thigh to ensure the correct size. washable and reusable

Code	Product	Pack size
3005	Ugo fix sleeve leg bag holder	Small 24-36cm
3006		Medium 36-55cm
3007		Large 40-70cm
3008		X Large 65-90cm
3027		XXLarge 75-105cm
10660S		Great Bear Libra leg bag holder
10660M	Medium 45-65cm	
10654L	Large 55-80 cm	
10665X	X Large 65-90cm	

Urinary catheter bed / 2litre bags Estimated usage 3 packs of 10 / monthly

Single Use 2 Litre bags (unsterile)

At night single use 2 litre bags should be attached to the bottom of the leg bag rather than drainable bags. The single use bags do not have any cost implications, when compared against drainable 2 litre bags **but are recommended by infection control** for use by all patients that have indwelling urinary catheters and require added drainage capacity overnight.

Code	Product	Pack size
10303V	Great Bear 2 litre non sterile night bag with non-return valve and single use T-tap slide action	10
PSU2	ProSys Non Sterile night drainage bag lever action	10

Drainable 2 Lt bags (sterile) Estimated usage 1 pack of 10 bags / 8 week

These 2 Litre catheter bags are for patients that are solely cared for in bed at all times, and do not use leg bag. They need to be changed every 5-7 days.

D813131	Bard Uriplan drainage bag, 2 litre with 98cm inlet tube, non-return valve and tap outlet.	10
DSBB2000	Bard Spirit Bed Bags, 2 litre with 120cm tube length, lever tap, 1 pair of latex free sterile gloves and needle free sample port	10

Night bag holders

Company	Tel no	These are not available on the NHS, but may be obtained free of charge from these companies if their products are being used.
Great Bear	0800 318559	
Bard	0800 4701930	
Coloplast	0800 374654	

Belly bags

Estimated usage 3 bags / 3 monthly

These are positioned on the abdomen. Urine drains into the bag by means of residual bladder pressure. Bag capacity is 1 litre. Bag needs to be changed every 28 days. Manufacturer can supply connector if night bag is required

There are other types of catheter bags available that are useful for patients with disabilities please speak to the Bladder and Bowel Service for advice before prescribing.

Code	Product	Pack size
B1000	Rusch Belly Bag 1000ml (short outlet tube)	1
B1000P	Rusch Belly Bag 1000ml (short outlet tube outlet with sample port)	1
B1000CT	Rusch Belly bag with extra long (60cm) outlet tube useful for wheelchair users	1

Miscellaneous supplies

Code	Product	Pack size	Estimated usage 10 straps as a one off prescription. Not to be added to repeat prescription
380812	Simpla leg bag straps	5 pairs	
10622C	Great Bear bag straps calf	10	
10637T	Great Bear bag straps leg	10	Estimated usage No more than 1 box every 6 months
	Sterile water for injection 10mls (plastic vials)	10	
MRB359	Clinipod	25 x 20mls	
62042	Conveen prep-wipes used prior to applying sheaths	54 wipes	
TR101	ConvaTec Niltac sting free medical adhesive remover	1X50mls	Needing to be prescribed as single items as indicated in pack size

10380A	Great Bear Soft Skin No Sting medical adhesive remover aerosol	1x50mls	
9811	Hollister Retracted Penis Pouch	10	

Intermittent Catheters

There are a wide range of intermittent catheters and therefore each patient needs to be assessed individually so that the correct catheter is supplied to meet their needs.

The Bladder and Bowel Service assess, initiate and make the ongoing supplies of catheters for patients intermittently self catheterising (ISC) . For those patients intermittently self catheterising but not presently under the care of the Bladder and Bowel Service, it is recommended that the GP contact the service in the first instance to discuss the review and ongoing management of patients ISC.

Please note: The prescription for intermittent catheters may need to be changed, either in quantity due to patient's needs, or as the patient becomes more skilled at intermittent self catheterisation the type of catheter requested may alter.

Urinary sheaths and bags

Prescribe 1 pack of 30 monthly

Sheaths are a useful alternative method to manage urinary incontinence for men.

- Sheaths must be accurately sized and the appropriate drainage system is prescribed.
- There are several different products listed (which are self-adhesive, silicone, latex free) to accommodate for good sizing both for girth fitting and length fitting.
- It is recommended that sheaths are changed every 24 hours to relieve pressure and to allow for personal hygiene.
- High quality breathable silicone for maximum comfort and skin protection.

Sheaths

Code	Product
Great Bear Libra Sheaths	
GBLSS24	Standard 24 mm
GBLSS28	Standard 28mm
GBLSS31	Standard 31 mm
GBLSS33	Standard 33 mm
GBLSS35	Standard 35 mm
GBLSS40	Standard 40mm
GBLSW24	Wide Band 24mm
GBLSW28	Wide Band 28mm
GBLSW31	Wide Band 31mm
GBLSW33	Wide Band 33 mm
GBLSW35	Wide Band 35 mm
GBLSW40	Wide Band 40mm

GBLSP24	Pop on 24 mm
GBLSP28	Pop on 28 mm
GBLSP31	Pop on 31 mm
GBLSP33	Pop on 33 mm
GBLSP35	Pop on 35 mm
GBLSP40	Pop on 40 mm
CliniSupplies Prosys Flofit Urinary Sheaths	
PSS24STD/B	Standard
PSS28STD/B	Standard
PSS31STD/B	Standard
PSS35STD/B	Standard
PSS40STD/B	Standard
PSS24WIDEBD/B	Wide band
PSS28WIDEBD/B	Wide Band
PSS31WIDEBD/D	Wide Band
PSS35WIDEBD/D	Wide Band
PSS40WIDEBD/D	Wide band
PSS24POPON/B	Pop on
PSS28POPON/B	Pop on
PSS31POPON/B	Pop on
PSS35POPON/B	Pop on
PSS40POPON/B	Pop on

Leg bag (Coloplast Conveen Security Plus) Prescribe 1 pack of 10 bags / 8 weekly
 Non-sterile adjustable anti-kink tubing leg bags for use with sheaths only.

Code	Product	Pack size
5164	350ml 30cm tubing non sterile	10
5165	350ml 50cm tubing non sterile	10
5160	500ml 30cm tubing non sterile	10
5161	500ml 50cm tubing non sterile	10

Small thigh bag (Conveen Active)

Useful for sheath use only when dribbling via the penis is experienced when a supra pubic is insitu or if dribbling incontinence is experienced some of the time and the patient doesn't want to use an incontinence pad all of the time.

Code	Product	Pack size
25501	250 ml Conveen Active 7cm adjustable inlet tube	10

Adhesive free, external collection system

For use by patients who experience some voiding control issues or are urinary incontinent, under advice of *Specialist Bladder and Bowel Nurse*
 Products issued with 1 instructional DVD

Afex system (iMEDicare):

Code	Product	Pack size
<i>A100-Plus brief code</i>	<i>Active Starter Pack- 1 High receptacle, 1 pair of briefs, 1 reusable collection bag (500ml 4-6 weeks), 1 cleaning container</i>	<i>1</i>
<i>A110-Plus brief code</i>	<i>Sedentary Starter Pack-1 low receptacle, 1 pair of briefs, 1 reusable collection bag (1200mls 4-6 weeks), 1 cleaning container, 1 cleansing solution</i>	<i>1</i>
<i>A120-Plus brief code</i>	<i>Mobility Assisted Starter Pack 1 supersoft (90 -120 days) receptacle, 1 pair of briefs, 1 reusable collection bag(1200mls 4-6 weeks), 1 cleaning container, 1 cleansing solution</i>	<i>1</i>

Afex briefs Prescribe 4 single units / yearly

Code	Product	Pack size
A200- S-08	Waist small 28"-30"	1
A200- M-08	Waist medium 32"-34	1
A200- L-08	Waist large 36"-38"	1
A200- XL-08	Waist XL 40"-42"	1
A200- 2X-08	Waist 2X 44"-46"	1
A200- 3X-08	Waist 3X 48"-50"	1

Afex urine collection bags Prescribe 9 bags / yearly

Code	Product	Pack size
A400-V-08	500ml vented bag	1
A400-B-08	500ml unvented bag	1
A400-L-08	1200ml leg bag	1
A400-E-08	Extra capacity	1
A400-X-08	2000ml bed bag	1

Afex receptacles Prescribe 3 units / yearly

Code	Product	Pack size
A300-L-08	Low receptacle	1
A300-S-08	Low supersoft Receptacle	1
A300-F-08	High Receptacle	1

Extension tubes and Stabiliser systems

Code	Product	Pack size
A700-A-08	Tube assembly complete (9 per year)	1
A700-S-07	Tube stabiliser	1
A710-S-09	Collection bag stabiliser-high grip	1

Anal Plugs (Coloplast)

Prescribe 1 pack of 20 / PRN

Patient needs to have an assessment by an appropriate healthcare professional prior to being prescribed this product. Anal plugs are used for patients with ano-rectal incontinence problems - most suitable for neurogenic conditions with impaired sensation. Plugs can be left in place up to 12 hours, if inserted correctly, with reasonable degree of comfort.

Code	Product	Pack size
1450	Peristeen Anal Plug Small	20
1451	Peristeen Anal Plug Large	20

Faecal collector (Hollister)

Prescribe 1 pack of 20 / PRN

Useful for bed bound patients experiencing semi-solid or liquid stools, protecting skin from irritating faecal discharge

Code	Product	Pack size
9822	50cm tubing non sterile 500ml	10
9821	50cm tubing non sterile 1000ml	10

Premier Pads for rectal bleeding

2004 (pack of 25) 2006 (pack of 15)	1X10cm X 20cm 1X20cm X 20cm	ESTIMATED USAGE Maximum quantity to be prescribed at any time one time = 100 pads
--	--------------------------------	---

Anal Irrigation

Irrigation Systems

Irrigation systems are cost effective when managed well and assessed by a continence advisor. It does stop hospital admissions, district nurses time.

If this option is to be considered for a patient, please refer the patient to either the Bladder and Bowel Service or Secondary Care for a specialist assessment prior to prescribing any equipment so that the patient is well informed about the procedure. If the patient is suitable for anal irrigation they will need to be taught the procedure which will normally be done in the community by the Bladder and Bowel Service. **Patients usually need support from the Bladder and Bowel Service when they initially use the system.**

There are four different types to meet the needs of patients.

Peristeen anal irrigation (Coloplast)

Code	Product	Pack size	<i>Please note:</i>
29121	<i>1 control unit, 2 rectal catheters, 1 water bag, 1 strap set, 1 toilet bag</i>	1	<i>Catheters are single use only.</i>
29126	<i>1 control unit, 2 small rectal catheters, 1 water bag, 1 strap set, 1 toilet bag</i>	1	<i>Water bag</i>

29122	15 Rectal catheters, 1 water bag	1	can only be used 15 times. Lid and control unit can only be used 90 times
29127	15 Small rectal catheters, 1 water bag	1	
29123	10 Rectal catheters	1	
29128	10 Small rectal catheters	1	
29124	1 Set of straps	1	
29125	2 Tubes and blue connectors	1	

Qufora IrriSedo System

Code	Product	Estimated Usage
<i>IrriSedo mini system 15 irrigation set 53601-015</i>	<i>IrriSedo Mini System (1 Pump, 15 hydrophilic coated cones 15 wast bags)</i>	<i>Prescribe 1 pack to last 1 month if wash out occurring on alternative days when established</i>
<i>IrriSedo Mini System 30 irrigation set 53601-030</i>	<i>IrriSedo Mini System 30 irrigation monthly set. (Pump, water bag, 15 cones, 15 waste bags and toilet bag)</i>	<i>Prescribe 1 pack to last 1 month if wash out occurring daily when established</i>
<i>IrriSedo Cone System 15 Irrigation Set QTM</i>	<i>(Contains 1 pump, 1 water bag, 15 hydrophilic coated cones, 15 waste bags)</i>	<i>Prescribe 1 pack to last 1 month if wash out occurring daily when established</i>
<i>IrriSedo Bed System (15 Irrigation Set) 53201-015</i>	<i>Contains 1 water bag with pump, 15 rectal catheters 12mm, 21mm 30mm 15 collection bags, 15 gel sachets)</i>	<i>Prescribe 1 pack to last 1 month if wash out occurring on alternative days when established</i>
<i>IrriSedo Balloon System - Base Set (Regular) 53101-015</i>	<i>(Contains 1 control unit with velcro strap, 1 water bag, 2 hydrophilic rectal balloon catheters (regular))</i>	<i>Prescribe 1 pack to last 1 month if wash out occurring on alternative days when established</i>

Braun IryPump Rectal Irrigation

A new system which involves instilling water into the bowel via the anus using a soft cone. The other benefit is the freedom just touching a button and this pumps the water into the rectum

29120E	<i>Irypump S Set</i>	<i>Quantity 1</i>
29220E	<i>Irycone (cone and tubing)</i>	<i>Accessories 1</i>
29240	<i>Irypump water container</i>	<i>Accessories 1</i>

NAVINA SMART

A new rectal irrigation system which would be useful for patients that have poor dexterity especially in spinal injury patients and neurogenic bladder and bowel

Code	Type	Content	Lifetime
69009	<i>Navina Smart system Regular for adult use only</i>	<i>1 navina smart control unit with power adaptor and cable. 1 water container, 1 Tube set, 2 regular rectal catheters, 1 accessory set, 1 navina case</i>	<i>400 uses or 2year whichever comes first control unit only</i>
69010	<i>Navian smart system small</i>	<i>1 navina smart control unit with power adaptor and cable. 1 water container, 1 Tube set, 2 small rectal catheters, 1 accessory set, 1 navina case</i>	<i>400 uses or 2year whichever comes first control unit only</i>
69003	<i>Consumable set regular for adult use only</i>	<i>1 water container, 15 regular Rectal catheters</i>	<i>15 uses Single use only</i>
69004	<i>Consumable set small</i>	<i>1 water container, 15 small Rectal catheters</i>	
68940	<i>Rectal catheter set Regular for adult use only</i>	<i>10 regular Rectal catheters</i>	
68943	<i>Rectal catheter set small</i>	<i>10 small Rectal catheters</i>	
69011	<i>Tube set</i>	<i>1 Catheter tube, 1 Water container tube</i>	<i>100 uses</i>
69014	<i>Navina smart control unit with power adaptor and cable</i>	<i>1 navina smart contrl unit 1 cable 400 uses or 2year whichever comes first control unit only 1 Power adaptor</i>	

Contenance appliances

Prescribe one only.

NB. Items with * need to have two litre drainable catheter bags also prescribed. There are a range of reusable appliances that can be prescribed that patients can use to enable them to remain continent even if their mobility is poor or they are chair or bed bound. Listed below are a range of products, but this is not extensive. Please contact the Bladder and Bowel Service for further advice. The contact details are at the front to this formulary.

Order code	Product	Pack size
6-50T	Beambridge male draining jug with tap*	1
6-26	Beambridge female saddle urinal	1
6-45	Beambridge Lady Jug *	1
6-18	Beambridge Bridge urinal	1
6-35	Beambridge Funnel (male) *	1
URIBag	Manfred Sauer Uribag (male)	1
URIBagF	Manfred Sauer Uribag (female)	1

Other useful information/links:

Bladder and Bowel Service: South Office (Worthing) 01273 265855

Bladder and Bowel Community: www.bladderandbowel.org

Association for Continence Advice www.aca.uk.com